

NEW! MRI Pads & Positioners



Cat# MFP-GNK

General MRI Kit Includes:

- two Rectangles 14 x 3 x 1
- two Rectangles 5 x 3 x 1
- two Rectangles 3 x 4 x 2
- two Wedges 10 x 10 x 3
- two Wedges 14 x 4 x 1.5
- two curve edge Wedges 5 x 5 x 3



Patients often complain about being uncomfortable during their MRI exam. A significant part of their discomfort relates to the mattress pads currently available for these MRI tables. Current table pads were not developed with the purpose of maximizing patient comfort nor do they utilize the latest technology that is available today. These New MRI Pads and Positioners are designed to conform to specific requirements for MRI compatibility and safety. Pads and positioners are compatible with existing tables and are designed to simply replace the current pads.

This unique visco-elastic material is specifically formulated to maximize patient comfort in a thin form factor to maximize free space inside MRI bores. This heat sensitive material molds to the patient's exact body shape and contours to evenly distribute weight and pressure. Features a Permanent Anti-microbial Agent with cover seams sown using coated nylon thread to minimize the risk of any biological contamination. Pad Covering is a Breathable Waterproof polymer fabric that is easy to maintain and clean with any standard disinfectant. Pads & Positioners are Resistant to stains and are Latex Free! Excellent radiolucent qualities. To increase safety, all pads have an anti-slip undersurface to keep the pads in place while patients are getting on and off the table.

MRI Kits Available for the following:

- | | |
|--------------|--|
| Cat# MFP-GE2 | MRI CTL Coil Pad Kit for GE 1.5T (2pc) |
| Cat# MFP-GE3 | MRI Table Pad Kit for GE 1.5T & 3.0T (3pc) |
| Cat# MFP-GE5 | MRI Table Pad Kit for GE 1.5T & 3.0T (5pc) |
| Cat# MFP-GE4 | MRI Table Pad Kit for GE .2T, .35, .7T (5pc) |
| Cat# MFP-GE7 | MRI Table Pad Kit for GE .2T, .35, .7T (3pc) |
| Cat# MFP-HT1 | MRI Table Pad Kit for Hitachi (8pc) |
| Cat# MFP-ST2 | MRI Table Pad Kit for Phillips, Seimens, & Toshiba (5pc) |
| Cat# MFP-PST | MRI Table Pad Kit for Phillips, Seimens, & Toshiba (3pc) |
| Cat# MFP-PT3 | MRI CTL Coil Pad Kit for Phillips, 1.5T (2pc) |

MRI Breast Biopsy Kit:

- | | |
|--------------|--|
| Cat# MFP-BBK | MRI Breast Biopsy Positioning Kit (10 pcs) |
|--------------|--|

OTHER POSITIONERS & PADS AVAILABLE.
CALL 1-800-251-2629 for pricing and information

NEW!

MRI Table Pads

Patients often complain about being uncomfortable during their MRI exam. A significant part of their discomfort relates to the mattress pads currently available for these MRI tables. Current table pads were not developed with the purpose of maximizing patient comfort nor do they utilize the latest technology that is available today. These New MRI Table Pads are designed to conform to specific requirements for MRI compatibility and safety. They are compatible with existing tables and are designed to simply replace the current pads.

This unique visco-elastic material is specifically formulated to maximize patient comfort in a thin form factor to maximize free space inside MRI bores. This heat sensitive material molds to the patient's exact body shape and contours to evenly distribute weight and pressure. Pads Feature a Permanent Antimicrobial Agent with cover seams sown using coated nylon thread to minimize the risk of any biological contamination. Pad Covering is a Breathable Waterproof polymer fabric that is easy to clean with standard disinfectant. Pads are Resistant to stains and are Latex Free! Excellent radiolucent quality. To increase safety, all pads have an anti-slip undersurface to keep the pads in place while patients are getting on/off the table.

Cat# MFP-20	MRI Table Pad for GE 1.5T & 3.0T 15"W x 13"L x 0.625"D
Cat# MFP-22	MRI Table Pad for GE 1.5T & 3.0T 15"W x 13"L x 1.25"D
Cat# MFP-24	MRI Table Pad for GE 1.5T & 3.0T 15"W x 27"L x 1.25"D (2pcs)
Cat# MFP-30	MRI Table Pad for GE .2T, .35, .7T 29"W x 14"L x 0.625"D
Cat# MFP-32	MRI Table Pad for GE .2T, .35, .7T 29"W x 14"L x 1.25"D
Cat# MFP-34	MRI Table Pad for GE .2T, .35, .7T 29"W x 28"L x 1.25"D (2pcs)
Cat# MFP-40	MRI Table Pad for Hitachi 31"W x 13.5"L x 0.625"D
Cat# MFP-42	MRI Table Pad for Hitachi 31"W x 27"L x 0.625"D (2pcs)
Cat# MFP-44	MRI Table Pad for Hitachi 27"W x 13.5"L x 0.625"D
Cat# MFP-46	MRI Table Pad for Hitachi 31"W x 13.5"L x 1.25"D
Cat# MFP-48	MRI Table Pad for Hitachi 31"W x 27"L x 1.25"D (2pcs)
Cat# MFP-50	MRI Table Pad for Phillips, Seimens, & Toshiba 19.5"W x 14"L x 0.625"D
Cat# MFP-52	MRI Table Pad for Phillips, Seimens, & Toshiba 19.5"W x 14"L x 1.25"D
Cat# MFP-54	MRI Table Pad for Phillips, Seimens, & Toshiba 19.5"W x 28"L x 1.25"D (2pcs)



**OTHER POSITIONERS & PADS AVAILABLE.
CALL 1-800-251-2629 for pricing and information**

