



**Cat# PLF**

## Radiation Attenuating Surgical Glove

Protection for your hands from direct and scatter radiation. For use during interventional procedures, diagnostic heart catheterization, coronary angioplasty, pain management, angiocardiology orthopedics and fluoroscopy procedures.

- Sterile - one time use
- 5 pair of one size per box
- Lead Free—no special disposal and non toxic for the environment
- Non Latex
- Powder Free to reduce the risk of allergic reaction
- High Tactile Sensitivity
- High Attenuation
- .18 mm Thickness
- Color: Black

### Determine your Glove Size:

*Glove sizing ranges from 6.5-9 in half sizes*

1. Measure around the widest point of the palm just below the knuckles of the dominant hand.
2. Match your measurement to the chart below for your estimated glove size.

### PROTECTION DATA:

60kVp	80kVp	100kVp	120kVp
56.5%	47.1%	41.6%	37.7%



*Test results in accordance with ASTM Test Method F2547-06*

SIZES	AVAILABLE
XS	6.5
S	7
M	7.5
L	8
XL	8.5
XXL	9

