

## **Techno-Aide**

## Axial/Sesamoid Podiatry Wedge Cat# YXNZ

High density 2 lb crosslink closed cell foam 15 degree wedge is an excellent positioning aide for viewing metatarsals, sesamoids, and plantar condyles. A smooth pleasant feel combined with superior physical and chemical properties makes cross-linked polyethylene foam a solution for an array of aplications that demand thicker sections of extremely fine cloed cell foam.

Surface is resistant to blood, Betadine, barium, alcohol, and other fluids. Antimicrobial: antibacterial and antifungal. Self deodorizing. Impervious to mildow, mold, and bacteria. Non toxic and contains no CFC's, HCFC's or hydrocarbon blowing agents.

Closedcell Foam size: 8 x 13 x 4



## WEIGHT BEARING DR/CR RECEPTOR HOLDER

Foam size: 17.5"L x 11.5"W x 4"H Slit to hold imaging receptor: .75" Wide

