

CT Mobile Shields

Bring a full 2.0mm Pb equivalency of radiation protection with our CT Protective Mobile Barriers

- Manufactured of the highest quality materials and craftsmanship
- 2.0mm Pb equivalency for in-room protection
- Transparent/clear upper section of x-ray protective Schott German-made shatterproof leaded glass Opaque lower section of solid steel with 2.0mm Pb Equiv. leaded x-ray protection.
- Locking caster system for easy mobility
- Perfect shielding for biopsy procedures.



Cat# PBMB-3060

Leaded Glass: 30" w x 60" h
 Opaque Panel: 30" w x 12" h
 Overall Size: 33" w x 75" h

Cat# PBMB-3024

Leaded Glass: 30" w x 24" h
 Opaque Panel: 30" w x 48" h
 Overall Size: 33" w x 75" h

Cat# PBMB-4830

Leaded Glass: 48" w x 30" h
 Opaque Panel: 48" w x 42" h
 Overall Size: 55" w x 75" h

*Crated for secure shipment. Shipping worldwide by ground, air or sea freight. Made in USA

