


NOTES:

FRAME MATERIAL:
STEEL WITH 2.0mm LEAD EQUIVALENCE


GLASS WINDOW:
2.0mm LEAD EQUIVALENCE GLASS

FRAME COLOR:
TAN


NOT TO SCALE



(A)
1 BARRIER ELEVATION



(B)
1 GLASS ELEVATION



(C)
1 SECTION AT FRAME



DWG No.	Rev.
PBMB-3024	
P/N:	
PBMB-3024	
DESC.	
Mobile Leaded Barrier w/24 inch Window	